Fabricant de la pomade :
Hersteller der Salbe:
Produttore della pomata:
Fabrikant van de zalf :
Manufacturer of the ointment:

Proizvođač masti:
Изготовитель мази
PIERRE FABRE MEDICAL DEVICES
29, avenue du Sidobre
81106 Castres Cedex
France

Fabricant des pansements :
Hersteller der Pflaster:
Produttore dei cerotti:
Fabrikant van de pleisters :
Manufacturer of the plasters:
Proizvođač flastera:
Изготовитель пластырей
FINESSE MEDICAL LTD
Unit 4, Royal Canal Business Park
Athlone Road
Longford
IRELAND

420647

FR

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EN

INSTRUCTION LEAFLET: Onyster®, ointment

What does Onyster® contain?

1 tube of 10g ointment containing:

Urea, lanolin, White soft paraffin, and 21 plasters

What is Onyster® used for?

Onyster® is recommended to be used as part of the topical treatment of fungal nail infections.

How does Onyster® work?

As a result of its composition that includes a 40% urea base, Onyster® weakens only affected areas of the nail (due to the fungus) responsible for the adhesion of the nail to the nail bed. This weakening allows selective detachment of the affected part of the nail, making it easy to cut and thus remove, in preparation for subsequent antifungal treatment.

How to use Onyster®?

- Apply the necessary amount of ointment to the infected nail once daily, covering the entire surface of the nail. Keep in place under an occlusive plaster for 24 hours.
- Repeat the application of the ointment and change the plaster every day, after removing the excess ointment from the

previous application. Depending on the extent of the infection and the thickness of the nail, 1 to 3 weeks of treatment may be necessary.

- After 1 to 3 weeks of treatment, it is recommended that a health care professional removes the detached, softened part of the nail using a pair of nail scissors.
- After the affected part of the nail has been removed, continue with an antifungal treatment

When should Onyster® not be used?

- If you are hypersensitive (allergic) to any of the components. In the absence of data, it is preferable not to use in combination with other products for topical use. Do not use without first consulting a health care professional.

Special instructions for use and warnings:

- In case of a hypersensitivity reaction (allergy), stop the application and inform your doctor.
- Do not swallow.
- Do not use on damaged skin.
- Wash your hands after use.
- Do not apply near the eyes or on mucous membranes.
- The dressings are for single use only, do not reuse.
- Do not use if the pack is damaged or if, on first opening, you notice that the tube does not have a protective seal.
- Keep out of the reach of children.

Possible side effects

This medical device may, in some individuals, cause side effects such as redness, irritation, softening, desquamation (flaking off of small skin cells from the skin surface), possibility of allergy, dryness of the skin.

Any other undesirable side effect should be reported to the manufacturer.

Special precaution for storage:

Store below 30°C.

CE mark year: 2008

Date of revision of the leaflet: March 2012

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