
Display description

7. Initial use

7.1 Inserting the batteries

1. Slide the battery compartment lid open.	2. Place the two batteries supplied in the pulse oximeter as shown. Ensure that the correct battery polarity is observed.	3. Close the battery compartment lid again.

7.2 Attaching the lanyard

To transport the pulse oximeter more easily you can attach a lanyard to the device.

1. Insert the narrow end of the lanyard through the holder as shown.	2. Draw the other end of the lanyard through the loop at the narrow end and tighten.

8. Operation

1. Insert one finger into the finger opening of the pulse oximeter as shown and hold it steady.	2. Press the function button. The pulse oximeter begins its measurement. Do not move during the measurement.	3. Your measurement values will appear on the screen after a few seconds.

- Note**
- When you remove your finger from the pulse oximeter, the device will automatically switch off after approx. five seconds.
- Function button**
- The function button on the pulse oximeter has two functions in total:
- Switch-on function:** When the pulse oximeter is switched off you can hold down the function button briefly to switch it on.
 - Brightness function:** To select your desired display brightness, hold down the function button for slightly longer during operation.

- Note:**
- The display orients automatically (vertical format, horizontal format). This ensures that the values are easy to read on the display at all times, regardless of how you hold the pulse oximeter.

9. Evaluating measurement results

WARNING	
The following table for evaluating your measurements does NOT apply to people with certain pre-existing conditions (e.g. asthma, heart failure, respiratory diseases) or whilst staying at altitudes above 1500 metres. If you have a pre-existing condition, always consult your doctor to evaluate your measurements.	
SpO ₂ measurement (oxygen saturation) in %	Classification/measures to be taken
99-94	Normal range
93-90	Decreased range: visit to the doctor recommended
< 90	Critical range Seek medical attention urgently

Source: based on Windisch W et al. S2k guideline: non-invasive and invasive ventilation as therapy for chronic respiratory failure Revision 2017; Pneumology 2017; 71: 722795

Decline in oxygen saturation depending on altitude

- Note**
- The following table informs you of the effects of various altitudes on oxygen saturation value and its impact on the human body. The following table does NOT apply to people with certain pre-existing conditions (e.g. asthma, heart failure, respiratory diseases etc.). People with pre-existing conditions can show signs of illness (e.g. hypoxia) at lower altitudes.

Altitude	Expected SpO ₂ value (oxygen saturation) in %	Impact on human body
1500-2500 m	> 90	No altitude sickness (normally)
2500-3500 m	~90	Altitude sickness, acclimatisation recommended
3500-5800 m	<90	Very frequent altitude sickness, acclimatisation absolutely essential
5800-7500 m	<80	Severe hypoxia, only limited length of stay possible
7500-8850 m	<70	Immediate, acute danger to life

Source: Hackett PH, Roach RC: High-Altitude Medicine. In: Auerbach PS (ed): Wilderness Medicine, 3rd edition; Mosby, St.Louis, MO 1995; 1-37.

10. Maintenance/cleaning

- IMPORTANT:**
- Do not use high-pressure sterilisation on the pulse oximeter! Under no circumstances should you hold the pulse oximeter under water, as this can cause liquid to enter and damage the pulse oximeter.
- Clean the housing and the interior rubber surface with a soft cloth dampened with medical alcohol after each use.
 - If a low battery status appears on the display of the pulse oximeter, change the batteries.
 - If you are not going to use the pulse oximeter for more than one month, remove both batteries from the device to avoid possible leaking.

11. Storage

- IMPORTANT:**
- Store the pulse oximeter in a dry place (relative humidity ≤95 %). If the humidity is too high it may shorten the service life of the pulse oximeter or damage it. Store the pulse oximeter in a place where the ambient temperature is between -40°C and 60°C.

- 12. Disposal**
- For environmental reasons, do not dispose of the device in the household waste at the end of its useful life. Dispose of the unit at a suitable local collection or recycling point. Dispose of the device in accordance with EC Directive – WEEE (Waste Electrical and Electronic Equipment). If you have any questions, please contact the local authorities responsible for waste disposal.
- The empty, completely flat batteries must be disposed of through specially designated collection boxes, recycling points or electronics retailers. You are legally required to dispose of the batteries.
- The codes below are printed on batteries containing harmful substances:
- Pb = Battery contains lead,
Cd = Battery contains cadmium,
Hg = Battery contains mercury.

13. What if there are problems?

Problem	Possible cause	Solution
The pulse oximeter is not displaying measurement values.	The batteries in the pulse oximeter are empty.	Replace the batteries.
	Batteries not inserted correctly.	Reinsert the batteries. If after reinserting the batteries correctly there are still no measurement values displayed, contact customer services.
The pulse oximeter is displaying measurement interruptions or high measurement value jumps.	Insufficient circulation in the measurement finger.	Observe the warnings and safety notes in chapter 5.
	Measurement finger is too large or too small.	Fingertip must have the following measurements: Width between 10 and 22 mm Thickness between 5 and 15 mm
	Finger, hand or body is moving.	Keep your finger, hand and body still during the measurement.
	Cardiac arrhythmia	Seek medical attention.

14. Technical Data

Type / Model.	Oxygen PO 01
Measurement method	Non-invasive measurement of arterial oxygen saturation of haemoglobin and pulse rate in finger
Measurement range	SpO ₂ 0 – 100%, Pulse 30 – 250 beats/minute
Accuracy	SpO ₂ 70 – 100%, ±2%, Pulse 30 – 250 bpm, ±2 beats/minute
Dimensions	L 61 mm x W 36 mm x H 32 mm
Weight	Approx. 58 g (including batteries)
Sensor to measure SpO ₂	Red light (wave length 660 nm); infra-red (wave length 880 nm); silicon receiver diode
Permissible operating conditions	+10°C to +40°C, <75 % relative humidity, 700–1060 hPa ambient pressure
Permissible storage conditions	-40°C to +60°C, ≤95 % relative humidity, 500–1060 hPa ambient pressure
Power supply	2x 1.5V — AAA batteries
Battery life	2 AAA batteries last for approx. 2 years of operation at 3 measurements per day (each of 60 seconds)
Classification	IP22, application part type BF

- The serial number is located on the device or in the battery compartment.
- Technical information is subject to change without notification to allow for updates.
- This device conforms with the European standards EN60601-1 and EN60601-1-2 (in accordance with CISPR 11, IEC 61000-4-2, IEC 61000-4-3, IEC 61000-4-8) and is subject to particular precautions with regard to electromagnetic compatibility. Please note that portable and mobile HF communication systems may interfere with this device. For more details, please contact our Customer Services at the address indicated.
 - This device complies with the EU Directive 93/42/EEC concerning medical devices, the Medizinproduktegesetz (German Medical Devices Act) and the DIN EN ISO 80601-2-61 standard (Medical electrical equipment – Particular requirements for the basic safety and essential performance of pulse oximeter equipment for medical use)

Notes on electromagnetic compatibility

- The device is suitable for use in all environments listed in these instructions for use, including domestic environments.
- The use of the device may be limited in the presence of electromagnetic disturbances. This could result in issues such as error messages or the failure of the display/device.
- Avoid using this device directly next to other devices or stacked on top of other devices, as this could lead to faulty operation. If, however, it is necessary to use the device in the manner stated, this device as well as the other devices must be monitored to ensure they are working properly.
- The use of accessories other than those specified or provided by the manufacturer of this device can lead to an increase in electromagnetic emissions or a decrease in the device's electromagnetic immunity; this can result in faulty operation.
- Keep portable RF communication devices (including peripheral equipment, such as antenna cables or external antennas) at least 30 cm away from all device parts, including all cables included in delivery. Failure to comply with the above can impair the performance of the device.
- Failure to comply with the above can impair the performance of the device.

15. Warranty/service

- Medel International, via Villapizzzone 26 / 20156, Milan (hereinafter referred to as “Medel”) provides a warranty for this product, subject to the requirements below and to the extent described as follows.
- The warranty conditions below shall not affect the seller's statutory warranty obligations which ensue from the sales agreement with the buyer. The warranty shall apply without prejudice to any mandatory statutory provisions on liability.**
- Medel guarantees the perfect functionality and completeness of this product.
- The worldwide warranty period is 5 years, commencing from the purchase of the new, unused product from the seller.
- The warranty only applies to products purchased by the buyer as a consumer and used exclusively for personal purposes in the context of domestic use.
- German law shall apply.
- During the warranty period, should this product prove to be incomplete or defective in functionality in accordance with the following provisions, Medel shall carry out a repair or a replacement delivery free of charge, in accordance with these warranty conditions.
- If the buyer wishes to make a warranty claim, they should approach their local retailer in the first instance: see the attached “International Service” list of service addresses.**
- The buyer will then receive further information about the processing of the warranty claim, e.g. where they can send the product and what documentation is required.
- A warranty claim shall only be considered if the buyer can provide Medel, or an authorised Medel partner, with
- a copy of the invoice/purchase receipt, and
 - the original product.
- The following are explicitly excluded from this warranty:
- deterioration due to normal use or consumption of the product;
 - accessories supplied with this product which are worn out or used up through proper use (e.g. batteries, rechargeable batteries, cuffs, seals, electrodes, light sources, attachments and nebuliser accessories);
 - products that are used, cleaned, stored or maintained improperly and/or contrary to the provisions of the instructions for use, as well as products that have been opened, repaired or modified by the buyer or by a service centre not authorised by Medel;
 - damage that arises during transport between manufacturer and customer, or between service centre and customer;
 - products purchased as seconds or as used goods;
 - consequential damage arising from a fault in this product (however, in this case, claims may exist arising from product liability or other compulsory statutory liability provisions).
- Repairs or an exchange in full do not extend the warranty period under any circumstances.
- Subject to errors and changes

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